



BD – Recall of PosiFlush™ SF Saline Flush Syringe

- On April 16, 2020, [BD announced](#) a voluntary, consumer-level medical device recall of multiple lots of [PosiFlush SF \(sterile field\)](#) saline flush because of holes in the packaging, which impacts package integrity and potentially compromises a sterile field.
 - While the sterility of the outer surface of the syringe may be compromised, the saline solution and the sterile path of the syringe are not impacted.
- BD distributed the affected recalled lots between February 28, 2019 through March 2, 2020.

Product Description	NDC#	Lot# (Expiration Date)
BD PosiFlush SF Saline Flush Syringe 10 mL	08290-3065-53	8353952 (12/31/21); 9011582 (12/31/21); 9017875 (12/31/21); 9024676 (01/31/22); 9045702 (01/31/22); 9060999 (02/28/22); 9079716 (02/28/22); 9127571 (02/28/22); 9143529 (04/30/22); 9156595 (05/31/22); 9163601 (05/31/22)

- PosiFlush SF saline flush is intended to be used only for the flushing of indwelling vascular access devices.
- When used in a sterile field, the compromised sterility due to holes in the packaging may increase the risk of infection to a patient, potentially leading to medical intervention and/or life-threatening injury. If the issue is identified prior to use and the syringe is discarded per standard clinical practice, this may lead to a delay or interruption of treatment and user dissatisfaction or annoyance.
- To date, BD has not received any adverse event reports related to this recall.
- If BD product is not available, customers may look for an equivalent sterile prefilled product in the marketplace.
- Lastly, as an alternate practice, clinicians may choose to draw up normal saline for a sterile field using a sterile syringe and sterile needle as described in the “Association of Surgical Technologists Guidelines for Safe Medication Practices in the Perioperative Area” and following their hospital guidelines.
- Anyone with an existing inventory of the recalled product should stop distribution and quarantine the product immediately.
- Contact BD by phone at **1-888-364-2985** for further information regarding this recall.



OptumRx® specializes in the delivery, clinical management and affordability of prescription medications and consumer health products. We are an Optum® company — a leading provider of integrated health services. Learn more at [optum.com](#).

All Optum® trademarks and logos are owned by Optum, Inc. All other brand or product names are trademarks or registered marks of their respective owners.

This document contains information that is considered proprietary to OptumRx and should not be reproduced without the express written consent of OptumRx.

RxNews® is published by the OptumRx Clinical Services Department.

©2020 Optum, Inc. All rights reserved.